INFRAEAZ

Infrared & Red Light Therapy Wrap

Swedish Posture®

EN The perfect combination wrap with 60 Leds using both Red Light and Infrared Light Therapy. Infrared Therapy (850Nm) is used for pain relief, muscle recovery and inflammation since it reaches deep into the cell tissues. Red Light Therapy (660Nm) targets the skin cells and helps with collagen production, wrinkles and aging spots.

Use InfraEaz when needed for 20 minute treatments. For optimal result use daily during 4 weeks. Use it directly on the skin.

NOTE that Infrared light is invisible to the human eye so it may appear that the LEDs are off when in fact the Infrared is on.

Led Qty: 60*0,5W

Rating Power: 30W

Battery duration: 115 min

Treatment time: 20 minutes

Battery charging time: 3 hours

USB-C charger

Wavelengths: 660Nm - 850Nm

Material: Neoprene

Low battery indication: Charging lamp starts flashing in red

Power button: Green light

- Press button once = Both Infrared and Red Light Therapy is on.
- Press button again = Only Red Light is on, Infrared off.
- Press button again = Only Infrared Light is on, Red light off (Infrared invisible).
- Press button again to turn off.

Power button: Blue light for pulsating therapy. Press button for more than 3 seconds.

- Press button once = Both Infrared and Red Light Therapy is on.
- Press button again = Only Red Light is on, Infrared off.
- Press button again = Only Infrared Light is on, Red light off (Infrared invisible)
- Press button again to turn off.

Use a damp towel to gently clean the wrap and elastic strap. Make sure the device is unplugged from any power and turned off. Only clean the fabric of the wrap and not the LED diodes directly. Do not submerge the wrap in water. Allow the wrap to air dry.

- Not intended for children and pregnant women.
- Do not exceed 30 minute treatments.
- Don't stare directly into the light.
- Never charge the belt unattended.

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Inspirera & Co AB/Swedish Posture